




We added to your treatment options by taking something away.





INTEGRA™ MATRIX WOUND DRESSING

Announcing Integra's New single-layer matrix wound dressing.



DESCRIPTION

INTEGRA[™] Matrix Wound Dressing is an advanced wound care device comprised of a porous matrix of cross-linked bovine tendon collagen and glycosaminoglycan. The collagen-glycosaminoglycan biodegradable matrix provides a scaffold for cellular invasion and capillary growth.

INDICATIONS

INTEGRA[™] Matrix Wound Dressing is indicated for the management of wounds including: partial and full-thickness wounds, pressure ulcers, venous ulcers, diabetic ulcers, chronic vascular ulcers, tunneled/undermined wounds, surgical wounds (donor sites/grafts, post-Moh's surgery, post-laser surgery, podiatric, wound dehiscence), trauma wounds (abrasions, lacerations, second-degree burns, skin tears) and draining wounds. The device is intended for one-time use.

HOW SUPPLIED

INTEGRA[™] Marix Wound Dressing is supplied sterile, in single use, double peel packages containing phosphate buffer. INTEGRA[™] Matrix Wound Dressing is available in the following sizes:

Product Codes	Sizes		Quantity
52021	2 inch x 2 inch	(5 cm x 5 cm)	1 unit/box
52025	2 inch x 2 inch	(5 cm x 5 cm)	5 units/box
54051	4 inch x 5 inch	(10 cm x 12.5 cm)	1 unit/box
54055	4 inch x 5 inch	(10 cm x 12.5 cm)	5 units/box
54101	4 inch x 10 inch	(10 cm x 25 cm)	1 unit/box
54105	4 inch x 10 inch	(10 cm x 25 cm)	5 units/box
58101	8 inch x 10 inch	(20 cm x 25 cm)	1 unit/box
58105	8 inch x 10 inch	(20 cm x 25 cm)	5 units/box

STORAGE

Store flat at room temperature. Avoid excessive heat (greater than 40° C). Avoid freezing.





Integra LifeSciences Corporation 311 Enterprise Drive • Plainsboro, NJ 08536 USA www.Integra-LS.com • www.ilstraining.com 1-877-444-1122 or 609-275-9004

INTEGRA Matrix Wound Dressing and the Integra wave logo are trademarks of Integra LifeSciences Corporation. © 2004 Integra LifeSciences Corporation. All Rights Reserved.